

# EuroVac

## LOW-PRESSURE POST-SURGERY WOUND DRAINAGE SYSTEM 100% CLOSED-CIRCUIT CONTINUOUS AND CONSTANT ASPIRATION

Revision date: 08.10.2025

- Aspiration pressure: 40 mmHg.
- Soft connector for inserting PVC and silicone drains from CH 10 to CH 18 (silicone drains require a special adaptor)
- Closed-circuit system with three non-return valves that make clamp handling unnecessary and promote hands-free operation.
- Extremely easy to use. Connect the drain and just compress the aspiration unit.
- Double hook for fixing the device to the patient's bed.
- The aspiration bag is emptied into the collection bag without interrupting the circuit.
- It is possible to remove the collection bag to allow the patient to deambulate more comfortably.
- Aspiration bag: volume 500 ml., graduated.
- Lower collection bag: volume 500 ml., graduated and equipped with an antibacterial filter.
- Aspiration bag: polyurethane  
Collection bag: PVC  
Non-return valves: silicone, PVC  
Spring: stainless steel  
Tubes: PVC, connectors: styrolux, caps: flexirene
- Latex free
- Double packaged (easy-peel pouch and inner bag)
- Sterilized by gamma rays
- Box of 25 pcs.



**CE**

**0425**

**CLASS Is  
RULE 1**

**STERILE R**



EuroVac's special non-return valves make clamp handling unnecessary, thus leaving the operator's hands free during the emptying procedure. Moreover, the absence of clamps eliminates the risk of clamp accidental opening leading to potentially hazardous fluid backflow in the system. So EuroVac is the safest and the most economical low-pressure drainage system on the market.

When the compressed aspiration unit has fully returned to its initial state, passive draining begins (when EuroVac is positioned below the wound level). In this way, the system's total capacity gains 500 ml of extra volume.

It is possible to remove the collection bag to allow the patient to deambulate more comfortably. The aspiration unit remains perfectly sealed, even with the lower bag removed.

**EV-S02** EuroVac 40mmHg, low-pressure closed-circuit disposable drainage with a removable collection bag.

**EV-S00** Removable collection bag.

## VACUUM CREATION / EMPTYING PROCEDURE

- Use both hands.
- Press the aspiration unit to the maximum level.
- During the emptying procedure the non-return valves prevent exudate backflow.



## COLLECTION BAG SUBSTITUTION

- To remove the bag twist it 90° anticlockwise.
- The new bag is positioned correctly when the arrows on the aspiration unit and the bag connectors are aligned.
- It is possible to remove the collection bag to save more space during deambulation.




Maximum capacity: 500 cc.

Initial aspiration capacity above wound level: 250 cc.

Initial aspiration capacity below wound level (passive mode): 350 cc. Tube length: 90 cm.

Aspiration bag size: 16 x 20 cm.

Collection bag diameter: 20 cm.



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy  
Tel. +39 051.817902 +39 051.817850  
[www.simitalia.it](http://www.simitalia.it) [www.medeurope.eu](http://www.medeurope.eu)  
E-Mail: [info@medeurope.eu](mailto:info@medeurope.eu)

Distributed by:

