


SIM Needle

Perforation needles



CLASS I

RULE 6



Description

SIM NEEDLE perforation needles are highly practical and convenient. Produced in AISI 304 steel, they can be sterilized and have an incorporated tube connector available in various diameters which permit a rapid and secure connection between the instrument and the drainage tube. The particular design of the cutting tip permits exceptional tissue perforation and considerably minimizes trauma.

NEEDLE CODE	NEEDLE CH	Needle for use with drainages having internal diameter not lower to:	Length needle mm.	MINIMUM QUANTITY pcs.
SN10160	9	2,5 mm.	160	5
SN12160	12	2,5 mm.	160	5
SN14160	15	3 mm.	160	5
SN16160	15	3,6 mm.	160	5
SN18160	18	3,7 mm.	160	5

SIM NEEDLE perforation needles are supplied in packages of 5 pcs. each.

WHILE STOCKS LAST

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
 Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845
 www.simitalia.it www.medeurope.eu
 E-Mail: info@medeurope.eu

Distributed by:

