

# SPRITZTUBE®

by MED EUROPE s.r.l.



[www.spritztube.eu](http://www.spritztube.eu)



[facebook.com/spritztube](https://facebook.com/spritztube)

## REVOLUTIONARY TRACHEAL CANNULA

STERILE DEVICE IN MEDICAL GRADE SILICONE



TWO SEPARATE INFLATING CUFFS, ONE FOR ESOPHAGUS OCCLUSION, ONE FOR SUPRAGLOTTIC POSITIONING

EASY TO USE

SECURE AND NON TRAUMATIC

VENTILATION AND INTUBATION WITH THE SAME DEVICE

CAN BE USED IN COMBINATION WITH ALL THE AVAILABLE INTUBATION METHODS (CLASSIC AND VIDEO LARYNGOSCOPE, BONFELS TYPE RIGID FIBERSCOPE AND FLEXIBLE FIBERSCOPE)



**EASY TO USE**

**FLEXIBLE**

**SECURE**



Via Gramsci, 9 - 40018 San Pietro in Casale (BO)

Tel.: +39.051.817902 - 051.810572 - Fax: +39.051.817845 - E-mail: [info@medeurope.eu](mailto:info@medeurope.eu) - Web: [www.medeurope.eu](http://www.medeurope.eu)

# SPRITZTUBE®

by MED EUROPE s.r.l.



[www.spritztube.eu](http://www.spritztube.eu)



[facebook.com/spritztube](https://facebook.com/spritztube)

## CANNULA TRACHEALE RIVOLUZIONARIA

PRESIDIO STERILE IN SILICONE MEDICAL GRADE


DUE DISTINTE CUFFIE DI GONFIAGGIO DEDICATE, UNA PER L'OCCLUSIONE ESOFAGEA ED UNA SOVRAGLOTTICA

FACILE DA UTILIZZARE

SICURO E NON TRAUMATICO

VENTILAZIONE ED INTUBAZIONE CON LO STESSO DISPOSITIVO

**PERMETTE L'INTUBAZIONE CON TUTTI I METODI DISPONIBILI**  
(LARINGOSCOPIO CLASSICO E A FIBRE OTTICHE, FIBROSCOPIO RIGIDO TIPO BONFILS, FIBROSCOPIO FLESSIBILE)



**FACILE DA USARE**

**FLESSIBILE**

**SICURO**



Via Gramsci, 9 - 40018 San Pietro in Casale (BO)

Tel.: +39.051.817902 - 051.810572 - Fax: +39.051.817845 - E-mail: [info@medeurope.eu](mailto:info@medeurope.eu) - Web: [www.medeurope.eu](http://www.medeurope.eu)