

Retractor for NEUROSURGERY

Double Hook for Galea

Revision date: 22.10.2020



Description


Double hook retractor for galea for use in neurosurgery.

Characteristics

- Made of AISI 301 steel
- Thickness 0,7 mm
- Length 33 mm
- "Medical grade" silicone elastic attachment
- Latex free device
- Packaged in double sterile envelope
- Sterilised by gamma rays
- Box of 50pz.

Code

REF MEDHOOK-005
(double hook)



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:

