




Endoperfect

The Ultimate System for Endoscopic Vessel Harvesting


ENDOSCOPIC VESSEL HARVESTING ADVANTAGES


- Fewer post-operative complications (fewer infections and less wound dehiscence)
- Less post-operative pain for the patient
- Shorter hospital stay
- Improved cosmetic results
- Reduction of the size of the surgical incision
- Uncompromised vessel quality
- No difference in operating time



Endoperfect

The Ultimate System
for Endoscopic Vessel Harvesting





EndoSafe

Safer for the patient

- Designed for maximum patient safety
- No need for closed insufflation: lower or no risk of CO2 embolism
- No risk of saphenous vein damage during vessel harvesting of the collateral vessels due to the special design of the manipulator.



EndoEasy

Easier to use

- The Endoperfect system is very easy and simple to use
- The Spacer device is designed to house the endoscopic instruments inside the operative channel which acts as a guide without the need for CO2 or manually lifting the instruments
- Optic housing channel separated from the spacer
- Smoke suction channel and lens washing channel to ensure the best possible visibility of the surgical field
- 30° optics: decreased risk of soiling when compared to a 0° optic



EndoLogic

Higher efficiency at a lower cost

- Reduced costs:
 - no need for specialized instrumentation
 - no need for expensive mono use optics or bipolar scalpels
- Versatile:
compatible with almost all currently available endoscopic devices



Endoperfect

technical
features

EndoPerfect

The Ultimate System for Endoscopic Vessel Harvesting



COD	DESCRIPTION	PACKAGING	PIECES PER PACKAGE
EP-001	EndoPerfect System for endoscopic vessel harvesting	n. 1 Tunnelling Device harvesting n. 1 Spacer/Dissector n. 1 Manipulator (Sterile Packs)	n. 1 EP-001 System

Patented product

Manufactured by:



Med. Europe
European Medical Supplies S.r.l.
Via Gramsci 9
40018 San Pietro in Casale (BO)
Tel. +39 051 810 572
Fax +39 051 817 850
e-mail: medeurope@simitalia.it

Distributed by:



Medical Instruments S.p.A.
Via Speranza 3/B
40068 San Lazzaro di Savena (BO)
Tel. +39 051 627 1111
Fax +39 051 627 1142
e-mail: info@medicalinstruments.it