

Med Stomach Tube

Gastric calibration tube



CE

0476

CLASS Is

RULE 5

STERILE R



The device is indicated for use in surgical procedures unique to the stomach. Introduced by mouth of the patient it allows, through its asymmetrical distal balloon, the movement of tissues, facilitating the surgeon involved in the intervention in laparoscopic or open surgery. The holes on the distal portion of atraumatic tip favor the complete suction.

Characteristics

- Material: medical grade silicone, shore 80.
- Length 800 mm.
- CH 38.
- Ext. diameter 12,5 mm
- Sterilised by gamma rays.
- Latex free.



REF MST-38

Packaging: 10 pieces/box

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
 Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845
 www.simitalia.it www.medeurope.eu
 E-Mail: info@medeurope.eu

Distributed by:

arfi