

SIM Silicon "Medical Grade"

Characteristics

Produced in pure "Medical Grade" silicone for multi-use medical/surgical purposes.

Water repellent. Maintains its flexibility over a wide range of temperatures and resists even closely repeated sterilizations. Shore 60.

Specifications

Non-sterile. Can be sterilized by autoclave at 121°C for 20 minutes or at 134°C for 6 minutes. Available in packages of 15 meter lengths in the following diameters:



REF	int. Ø mm.	ext. Ø mm.	CH
1x3	1	3	9
1,02x2,16	1,02	2,16	7
1,57x2,41	1,57	2,41	7
1,57x3,18	1,57	3,18	10
1,98x3,18	1,98	3,18	10
2,64x4,88	2,64	4,88	13
3x4	3	4	12
3,18x6,35	3,18	6,35	19
3,35x4,65	3,35	4,65	14
4,8x7,9	4,8	7,9	24
4,8x9,5	4,8	9,5	29

REF	int. Ø mm.	ext. Ø mm.	CH
5x7	5	7	21
5,8x11,18	5,8	11,18	34
6x8	6	8	24
6,4x9,5	6,4	9,5	29
6,4x11,1	6,4	11,1	33
6,4x12,7	6,4	12,7	38
6,9x9,5	6,9	9,5	28
7x12	7	12	36
7,9x12,7	7,9	12,7	38
8x12	8	12	36
9,5x14,3	9,5	14,3	43



0476

CLASS IIa

RULE 7



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
 Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845
 www.simitalia.it www.medeurope.eu
 E-Mail: info@medeurope.eu

Distributed by:

