




Aspiration Kit

Steel suction cannula complete with drop shaped suction control handle and connection tube

- By its drop shaped suction interrupter, FUKUSHIMA -Design suction cannula ensures excellent fine control of the suction power
- Rounded edges, atraumatic distal end of lumen
- Cannulae are made of Steel AISI 304
- Complete with silicone tube, length 60cm., equipped with rigid plastic fixation cone at one end and silicone funnel at the other end.
- Internal diameter is uniform over the entire length of the cannula
- Suction control handle made of SBC plastic material
- Matches all commonly used connection tubes
- Colour Code for prompt recognition
- Latex free
- Sterilised by gamma radiation
- Double sterile peel packaging



REF Cannula + Connection tube	Internal Diameter mm.	External Diameter mm.	CH	Cannula Length (mm)	Handle colour	Pcs. Per BOX
KIT G01*	0,96	1,4	4	90	WHITE	20
KIT G02*	1,45	2,0	6	80	GREEN	20
KIT G03*	2,54	3,0	9	80	BLUE	20
KIT G04*	3,0	4,0	12	80	WHITE	20
KIT G05*	0,96	1,4	4	130	WHITE	240
KIT G06	1,45	2,0	6	130	GREEN	20
KIT G07	2,54	3,0	9	130	BLUE	20
KIT G08	3,0	4,0	12	130	WHITE	20
KIT G09	1,45	2,0	6	160	GREEN	20
KIT G10	2,54	3,0	9	160	BLUE	20
KIT G11	3,0	4,0	12	160	WHITE	20

*Product not in stock, available by specific request

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu




Distributed by:



Aspiration Kit

Steel suction cannula complete with drop shaped suction control handle and connection tube

- There are 6 distal lateral holes to reduce tissue trauma
- By its drop shaped suction interrupter, FUKUSHIMA-Design suction cannula ensures excellent fine control of the suction power
- Rounded edges, atraumatic distal end of lumen
- Cannulae are made of Steel AISI 304
- The kit is completed with silicone tube, length 60cm., equipped with rigid plastic fixation cone at one end and silicone funnel at the other end
- Internal diameter is uniform over the entire length of the cannula
- Suction control handle made of SBC plastic material
- Matches all commonly used connection tubes
- Colour Code for prompt recognition
- Latex free
- Sterilised by gamma radiation
- Packaging: outer sterile blister and inner protective bag (single cannula); double sterile pouch (kit).



REF Cannula	Internal Diameter (mm)	External Diameter (mm)	CH	Cannula Length (mm)	Handle colour	Pcs. per BOX
CSSR20F150*	1,45	2,0	6	150	GREEN	60
CSSR30F150*	2,54	3,0	9	150	BLUE	60
Cannula + Connection tube	Internal Diameter (mm)	External Diameter (mm)	CH	Cannula Length (mm)	Handle colour	Pcs. per BOX
KIT G12*	0,96	1,4	4	150	WHITE	20
KIT G13*	1,45	2,0	6	150	GREEN	20
KIT G14*	2,54	3,0	9	150	BLUE	20

* Product not in stock, available by specific request

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:

