



Kit for Middle Ear Surgery

Revision date: 18.04.2023



- The kit comes ready for use and cannulae are opportunely placed to avoid damaging the package
- Non-traumatic rounded tips
- Cannulae made of AISI 304 medical grade steel
- Uniform internal diameter over the entire cannula length
- Angled
- Suction control handle made of SBC plastic material
- PVC connection tube
- Colour coded for quick recognition
- Latex free
- Sterilised by gamma rays
- Packaged in single sterile envelope



Kit includes:

- 1 cannula 90mm. length, external diameter 1.0mms., CH3
- 1 cannula 90mm. length, external diameter 1.4mms., CH4
- 1 cannula 80mm. length, external diameter 2.0mms., CH6
- 1 connection tube (PVC), length 3500mms., external diameter 8.0 mms.

Complete kit REF: KIT002*

Box of 40 kits (4 boxes of 10 kits each)

*Product not in stock. Made to order upon request.

Manufactured by:




Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850 www.simitalia.it
www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:



Kit for Stapedectomy and Middle Ear Surgery

Revision date: 18.04.2023



- The kit comes ready for use and cannulae are opportunely placed to avoid damaging the package
- Non-traumatic rounded tips
- Cannulae made of AISI 304 medical grade steel
- Uniform internal diameter over the entire cannula length
- Angled
- Suction control handle made of SBC plastic material
- PVC connection tube
- Colour coded for quick recognition
- Latex free
- Sterilised by gamma rays
- Packaged in single sterile envelope



Kit includes:

- 2 cannulae length 90mms., external diameter 0.7mms., CH2
- 2 cannulae length 90mms., external diameter 1.0mms., CH3
- 2 cannulae length 90mms., external diameter 1.4mms., CH4
- 2 cannulae length 80mms., external diameter 2.0mms., CH6
- 1 connection tube (PVC), length 3500mms., external diameter 8.0 mms.

Complete kit REF: KIT001*

Box of 40 kits (4 boxes of 10 kits each)

*Product not in stock. Made to order upon request.

Manufactured by:




Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:




Kit for Middle Ear Surgery

Revision date: 18.04.2023



- The kit comes ready for use and cannulae are opportunely placed to avoid damaging the package
- Non-traumatic rounded tips
- Cannulae made of AISI 304 medical grade steel
- Uniform internal diameter over the entire cannula length
- Angled
- Suction control handle made of SBC plastic material
- PVC connection tube
- Colour coded for quick recognition
- Latex free
- Sterilised by gamma rays
- Packaged in single sterile envelope



Kit includes:

- 1 cannulae length 90mm., external diameter 0.7mm., CH2
- 1 cannulae length 90mm., external diameter 1.0mm., CH3
- 1 cannulae length 80mm., external diameter 4.0mm., CH12
- 1 connection tube (PVC), length 3500mms., external diameter 8.0 mms.

REF KIT004* Box of 40 kits (4 boxes of 10 kits each)

*Product not in stock. Made to order upon request.

Manufactured by:




Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:



Kit for Middle Ear Surgery

Revision date: 18.04.2023



- The kit comes ready for use and cannulae are opportunely placed to avoid damaging the package
- Non-traumatic rounded tips
- Cannulae made of AISI 304 medical grade steel
- Uniform internal diameter over the entire cannula length
- Angled
- Suction control handle made of SBC plastic material
- PVC connection tube
- Colour coded for quick recognition
- Latex free
- Sterilised by gamma rays
- Packaged in single sterile envelope



Kit includes:

- 1 cannulae length 90mm., external diameter 0.7mm., CH2
- 1 cannulae length 90mm., external diameter 1.0mm., CH3
- 1 cannulae length 90mm., external diameter 1.4mm., CH4
- 1 cannulae length 80mm., external diameter 2.0mm., CH6
- 1 cannulae length 80mm., external diameter 3.0mm., CH9
- 1 connection tube (PVC), length 3500mms., external diameter 8.0 mms.

REF KIT003* Box of 40 kits (4 boxes of 10 kits each)

*Product not in stock. Made to order upon request.

Manufactured by:




Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:



Kit for Middle Ear Surgery

Revision date: 23.04.2026



- The kit comes ready for use and cannulae are opportunely placed to avoid damaging the package
- Non-traumatic rounded tips
- Cannulae made of AISI 304 medical grade steel
- Uniform internal diameter over the entire cannula length
- Angled
- Suction control handle made of SBC plastic material
- PVC connection tube
- Colour coded for quick recognition
- Latex free
- Sterilised by gamma rays
- Packaged in single sterile envelope



Kit includes:

- 1 cannulae length 160mm., external diameter 2.0mm., CH6
- 1 cannulae length 160mm., external diameter 3.0mm., CH9
- 1 cannulae length 160mm., external diameter 3.0mm., CH9
- 1 cannulae length 130mm., external diameter 4.0mm., CH12
- 1 cannulae length 140mm., external diameter 4.0mm., CH12
- 1 cannulae length 160mm., external diameter 4.0mm., CH12
- 1 connection tube (PVC), length 3500mms., external diameter 8.0 mms.

REF KIT006* Box of 40 kits (4 boxes of 10 kits each)

*Product not in stock. Made to order upon request.

Manufactured by:




Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:



Kit for Middle Ear Surgery

Revision date: 23.04.2026



- The kit comes ready for use and cannulae are opportunely placed to avoid damaging the package
- Non-traumatic rounded tips
- Cannulae made of AISI 304 medical grade steel
- Uniform internal diameter over the entire cannula length
- Angled
- Suction control handle made of SBC plastic material
- PVC connection tube
- Colour coded for quick recognition
- Latex free
- Sterilised by gamma rays
- Packaged in single sterile envelope



Kit includes:

- 1 cannulae length 160mm., external diameter 3.0mm., CH9
- 1 cannulae length 160mm., external diameter 3.0mm., CH9
- 1 cannulae length 140mm., external diameter 4.0mm., CH12
- 1 cannulae length 160mm., external diameter 4.0mm., CH12
- 1 connection tube (PVC), length 3500mms., external diameter 8.0 mms.

REF KIT007* Box of 40 kits (4 boxes of 10 kits each)

*Product not in stock. Made to order upon request.

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:

