


MEDRETRACTOR—MEDHOOK SURGICAL TISSUE RETRACTOR

Steel hook for traction or tissue dilation during surgery

Revision date: 02.07.2026



Description

A retractor consisting of a steel hook complete with a silicone elastic band for fixing.

Used in surgery to separate tissue flaps following incisions and allow greater exposure of the underlying organs in the surgical field.

Characteristics

- Made of AISI 316 steel.
- Diameter 1,6mm.
- “Medical grade” silicone elastic band for fixing.
- Latex free device.
- Packaged in double sterile pouch.
- Gamma ray sterilization.
- Box of 50 pcs.



Codes

REF **MEDHOOK001**

Large size, 110mms.

REF **MEDHOOK002**


Large size, 110mms., coated with “medical grade” silicone.

REF **MEDHOOK003**

Small size, 100mms.

REF **MEDHOOK004**

Small size, 100mms., coated with “medical grade” silicone



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:



MEDRETRACTOR—MEDHOOK SURGICAL TISSUE RETRACTOR

Double steel hook for traction or tissue dilation during surgery


Revision date: 02.07.2026



Description

A double retractor consisting of a steel hook complete with a silicone elastic band for fixing.

Used in surgery to separate tissue flaps following incisions and allow greater exposure of the underlying organs in the surgical field.




Characteristics

- Made of AISI 301 steel.
- Diameter 0,7 mm.
- Length 33 mm.
- "Medical grade" silicone elastic band for fixing.
- Latex free device.
- Packaged in double sterile pouch.
- Sterilised by gamma rays.
- Box of 50 pcs.

Code

REF **MEDHOOK005**

Double hook



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:



MEDRETRACTOR—GIMIHOOO

SURGICAL TISSUE RETRACTOR

Steel hook for traction or tissue dilation during surgery

Revision date: 31.08.2026



CLASS IIa

RULE 7

STERILE R



Description

A retractor consisting of a steel hook complete with a silicone elastic band for fixing. Used in surgery to separate tissue flaps following incisions and allow greater exposure of the underlying organs in the surgical field.

Characteristics

- Made of AISI 316 steel.
- “Medical grade” silicone elastic band with different hardness (20 ShA, 40 ShA, 60 ShA).
- Available in large or small versions, with parallel or twisted wire.
- Latex free device.
- Sterilised by gamma rays.
- Double sterile packaging: inner blister and outer pouch.
- Box of 20 pcs.

Codes

REF UNCS20P Small hook with 20ShA silicone elastic band
REF UNCS40P Small hook with 40ShA silicone elastic band
REF UNCS60P Small hook with 60ShA silicone elastic band
REF UNCS20G Large hook with 20ShA silicone elastic band
REF UNCS40G Large hook with 40ShA silicone elastic band
REF UNCS60G Large hook with 60ShA silicone elastic band
REF UNCR20P Small twisted hook with ShA 20 silicone elastic band
REF UNCR40P Small twisted hook with ShA 40 silicone elastic band
REF UNCR60P Small twisted hook with ShA 60 silicone elastic band
REF UNCR20G Large twisted hook with ShA 20 silicone elastic band
REF UNCR40G Large twisted hook with ShA 40 silicone elastic band
REF UNCR60G Large twisted hook with ShA 60 silicone elastic band



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simalia.it www.medeurope.eu

E-Mail: info@medeurope.eu


Distributed by:



MEDRETRACTOR—GIMIHOO

SURGICAL TISSUE RETRACTOR

Set of steel hooks for traction or tissue dilation during surgery



Description

A retractor consisting of a set of two steel hook complete with silicone elastic band for fixing and a polyamide clamp for anchoring the retractor system to the operating table.

Used in surgery to separate tissue flaps following incisions and allow greater exposure of the underlying organs in the surgical field.

Characteristics

- Made of AISI 316 steel.
- “Medical grade” silicone elastic band and polyamide fixing clamp.
- Latex free device.
- Double sterile packaging: inner blister and outer pouch.
- Sterilised by gamma rays.
- Box of 20 pcs.



CLASS IIa
RULE 7



Codes

REF SETUNCS20P Set of two small hooks with 20 ShA silicone elastic band and fixing clamp
 REF SETUNCS40P Set of two small hooks with 40 ShA silicone elastic band and fixing clamp
 REF SETUNCS60P Set of two small hooks with 60 ShA silicone elastic band and fixing clamp
 REF SETUNCS20G Set of two large hooks with 20 ShA silicone elastic band and fixing clamp
 REF SETUNCS40G Set of two large hooks with 40 ShA silicone elastic band and fixing clamp
 REF SETUNCS60G Set of two large hooks with 60 ShA silicone elastic band and fixing clamp
 REF SETUNCR20P Set of two small twisted hooks with 20 ShA silicone elastic band and fixing clamp
 REF SETUNCR40P Set of two small twisted hooks with 40 ShA silicone elastic band and fixing clamp
 REF SETUNCR60P Set of two small twisted hooks with 60 ShA silicone elastic band and fixing clamp
 REF SETUNCR20G Set of two large twisted hooks with 20 ShA silicone elastic band and fixing clamp
 REF SETUNCR40G Set of two large twisted hooks with 40 ShA silicone elastic band and fixing clamp
 REF SETUNCR60G Set of two large twisted hooks with 60 ShA silicone elastic band and fixing clamp

Accessories Codes

UNCLAMP Clamp for hooks
 RET20 20 ShA silicone elastic band for retractor
 RET40 40 ShA silicone elastic band for retractor
 RET60 60 ShA silicone elastic band for retractor

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:

