


MEDRETRACTOR—MEDHOOK SURGICAL TISSUE RETRACTOR

Steel hook for traction or tissue dilation during surgery

Revision date: 02.07.2026



Description

A retractor consisting of a steel hook complete with a silicone elastic band for fixing. Used in surgery to separate tissue flaps following incisions and allow greater exposure of the underlying organs in the surgical field.

Characteristics

- Made of AISI 316 steel.
- Diameter 1,6mm.
- “Medical grade” silicone elastic band for fixing.
- Latex free device.
- Packaged in double sterile pouches.
- Gamma ray sterilization.
- Box of 50 pcs.



Codes

REF MEDHOOK001

Large size, 110mms.

REF MEDHOOK002


Large size, 110mms., coated with “medical grade” silicone.

REF MEDHOOK003

Small size, 100mms.

REF MEDHOOK004

Small size, 100mms., coated with “medical grade” silicone.



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:



MEDRETRACTOR—MEDHOOK SURGICAL TISSUE RETRACTOR

Double steel hook for traction or tissue dilation during surgery

Revision date: 02.07.2026



Description

A double retractor consisting of a steel hook complete with a silicone elastic band for fixing.

Used in surgery to separate tissue flaps following incisions and allow greater exposure of the underlying organs in the surgical field.



CLASS IIa

RULE 7

0425

STERILE R




Characteristics

- Made of AISI 301 steel.
- Diameter 0,7 mm.
- Length 33 mm.
- "Medical grade" silicone elastic band for fixing.
- Latex free device.
- Packaged in double sterile pouches.
- Sterilised by gamma rays.
- Box of 50 pcs.

Code

REF MEDHOOK005

Double hook



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu

E-Mail: info@medeurope.eu


Distributed by:



MEDRETRACTOR—GIMIHOO

SURGICAL TISSUE RETRACTOR

Steel hook for traction or tissue dilation during surgery



CLASS IIa

RULE 7

STERILE R



Description

A retractor consisting of a steel hook complete with a silicone elastic band for fixing.

Used in surgery to separate tissue flaps following incisions and allow greater exposure of the underlying organs in the surgical field.

Characteristics

- Made of AISI 316 steel.
- “Medical grade” silicone elastic band with different hardness (20 ShA, 40 ShA, 60 ShA).
- Available in large or small versions, with parallel or twisted wire.
- Latex free device.
- Packaged in internal blister and external pouch.
- Sterilised by gamma rays.
- Box of 20 pcs.



Codes

REF UNCS20P Small hook with 20ShA silicone elastic band

REF UNCS40P Small hook with 40ShA silicone elastic band

REF UNCS60P Small hook with 60ShA silicone elastic band

REF UNCS20G Large hook with 20ShA silicone elastic band

REF UNCS40G Large hook with 40ShA silicone elastic band

REF UNCS60G Large hook with 60ShA silicone elastic band

REF UNCR20P Small twisted hook with ShA 20 silicone elastic band

REF UNCR40P Small twisted hook with ShA 40 silicone elastic band

REF UNCR60P Small twisted hook with ShA 60 silicone elastic band

REF UNCR20G Large twisted hook with ShA 20 silicone elastic band

REF UNCR40G Large twisted hook with ShA 40 silicone elastic band

REF UNCR60G Large twisted hook with ShA 60 silicone elastic band



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu

E-Mail: info@medeurope.eu


Distributed by:



MEDRETRACTOR—GIMIHOOO

SURGICAL TISSUE RETRACTOR

Set of steel hook for traction or tissue dilation during surgery




Description

A retractor consisting of a set of two steel hook complete with silicone elastic band for fixing and a polyamide clamp for anchoring the retractor system to the operating table.

Used in surgery to separate tissue flaps following incisions and allow greater exposure of the underlying organs in the surgical field.

Characteristics

- Made of AISI 316 steel.
- “Medical grade” silicone elastic band and polyamide fixing clamp.
- Latex free device.
- Packaged in internal blister and external pouch.
- Sterilised by gamma rays.
- Box of 20 pcs.



Codes

REF SETUNCS20P Set of two small hooks with 20 ShA silicone elastic band and fixing clamp

REF SETUNCS40P Set of two small hooks with 40 ShA silicone elastic band and fixing clamp

REF SETUNCS60P Set of two small hooks with 60 ShA silicone elastic band and fixing clamp

REF SETUNCS20G Set of two large hooks with 20 ShA silicone elastic band and fixing clamp

REF SETUNCS40G Set of two large hooks with 40 ShA silicone elastic band and fixing clamp

REF SETUNCS60G Set of two large hooks with 60 ShA silicone elastic band and fixing clamp

REF SETUNCR20P Set of two small twisted hooks with 20 ShA silicone elastic band and fixing clamp

REF SETUNCR40P Set of two small twisted hooks with 40 ShA silicone elastic band and fixing clamp

REF SETUNCR60P Set of two small twisted hooks with 60 ShA silicone elastic band and fixing clamp

REF SETUNCR20G Set of two large twisted hooks with 20 ShA silicone elastic band and fixing clamp

REF SETUNCR40G Set of two large twisted hooks with 40 ShA silicone elastic band and fixing clamp

REF SETUNCR60G Set of two large twisted hooks with 60 ShA silicone elastic band and fixing clamp

Accessories Codes

UNCLAMP Clamp for hooks

RET20 20 ShA silicone elastic band for retractor

RET40 40 ShA silicone elastic band for retractor

RET60 60 ShA silicone elastic band for retractor

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:

