


SIM MICRO CANNULA “Slow Reaction”

Steel suction cannula with suction control handle



Why?

- Because reusable suction cannulae with fine lumens are very difficult to clean
- Because disposable suction cannulae are recommended for very high risk procedures, especially neuro and ophthalmic.

Characteristics

- Aspiration control hole has a diameter of 1,5 mm., allowing the surgeon a slower aspiration response. Whenever the hole is opened the aspiration on the tip will slowly decrease. Even with a very high aspiration pressure the cannula makes only a low amount of noise.
- Complete line and specials can be easily produced following your indications
- Rounded edges, atraumatic distal end of lumen
- Non blocking design
- Cannulae are made of Steel AISI 304
- Single angle
- Internal diameter is uniform over the entire length of the cannula
- Suction control handle made of SBC plastic material.
- Matches all commonly used connection tubes
- Colour code for prompt recognition. Latex free
- Sterilised by gamma radiation. Single sterile blister packaging.



0425

CLASS IIa

RULE 6

STERILE R



Code	Internal Diameter mm.	External Diameter mm.	CH	Cannula length mm.	Handle colour	Pcs. per BOX	Applications
CSSR07X70	0,43	0,7	2	70	BLACK	60	ENT
CSSR07X90	0,43	0,7	2	90	BLACK	60	ENT
CSSR10X70	0,7	1	3	70	BROWN	60	ENT
CSSR10X90	0,7	1	3	90	BROWN	60	ENT
CSSR14X70	0,96	1,4	4	70	WHITE	60	ENT
CSSR14X90	0,96	1,4	4	90	WHITE	60	ENT
CSSR20X80	1,45	2,0	6	80	GREEN	60	ENT, Plastics, Otological Surgery & outpatients
CSSR20X130	1,45	2,0	6	130	GREEN	60	General, Plastics, Obs & Gynae, Paediatrics, Neuro, Oral & Max Fax
CSSR20X160	1,45	2,0	6	160	GREEN	60	General, Plastics, Obs & Gynae, Paediatrics, Neuro, Oral & Max Fax
CSSR30X80	2,54	3,0	9	80	BLUE	60	ENT, Plastics
CSSR30X130	2,54	3,0	9	130	BLUE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics
CSSR30X160	2,54	3,0	9	160	BLUE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics
CSSR40X80	3,0	4,0	12	80	WHITE	60	ENT, Plastics
CSSR40X130	3,0	4,0	12	130	WHITE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics
CSSR40X160	3,0	4,0	12	160	WHITE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu


E-Mail: info@medeurope.eu

Distributed by:



SIM MICRO CANNULA “Slow Reaction”

Steel suction cannula with suction control handle and plastic cleaning stylet



0425

CLASS IIa

RULE 6

STERILE R



Why?

- Because reusable suction cannulae with fine lumens are very difficult to clean
- Because disposable suction cannulae are recommended for very high risk procedures, especially neuro and ophthalmic

Characteristics

- Aspiration control hole has a diameter of 1,5 mm., allowing the surgeon a slower aspiration response. Whenever the hole is opened the aspiration on the tip will slowly decrease. Even with a very high aspiration pressure the cannula makes only a low amount of noise
- Equipped with a stylet in plastic material for cleaning the cannula in case of obstruction
- Complete line and specials can be easily produced following your indications
- Rounded edges, atraumatic distal end of lumen
- Non blocking design
- Cannulae are made of Steel AISI 304
- Single angle
- Internal diameter is uniform over the entire length of the cannula
- Suction control handle made of SBC plastic material
- Matches all commonly used connection tubes
- Colour code for prompt recognition
- Latex free
- Sterilised by gamma radiation
- Single sterile blister packaging



* Product not in stock, available by specific request

Code	Internal Diameter mm.	External Diameter mm.	CH	Cannula length mm.	Handle colour	Pcs. per BOX	Applications
CSSR30S80*	2,54	3,0	9	80	BLUE	60	ENT, Plastics
CSSR30S130*	2,54	3,0	9	130	BLUE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics
CSSR30S140*	2,54	3,0	9	140	BLUE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics
CSSR30S160*	2,54	3,0	9	160	BLUE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics
CSSR40S140*	3,0	4,0	12	140	WHITE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics
CSSR40S160*	3,0	4,0	12	160	WHITE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:



SIM MICRO CANNULA “Slow Reaction” 90°L

Steel suction cannula rotated 90° to the left

Complete with suction control handle

Why?

- Because reusable suction cannulae with fine lumens are very difficult to clean
- Because disposable suction cannulae are recommended for very high risk procedures, especially neuro and ophthalmic

Characteristics

- Aspiration control hole has a diameter of 1,5 mm., allowing the surgeon a slower aspiration response. Whenever the hole is opened the aspiration on the tip will slowly decrease. Even with a very high aspiration pressure the cannula makes only a low amount of noise.
- Complete line and specials can be easily produced following your indications
- Rounded edges, atraumatic distal end of lumen
- Non blocking design
- Cannulae are made of Steel AISI 304
- Single angle
- Internal diameter is uniform over the entire length of the cannula
- Angle and rotate to 90° in relation to the handle to facilitate a lateral hold
- Suction control handle made of SBC plastic material
- Matches all commonly used connection tubes
- Colour Code for prompt recognition
- Latex free
- Sterilised by gamma radiation
- Single sterile blister packaging





CLASS IIa

0425

RULE 6

STERILE R



Code	Int. Diameter mm.	Ext. Diameter mm.	CH	Cannula length mm.	Handle colour	Pcs. per BOX
CSSR20L130	1,45	2,0	6	130	GREEN	60
CSSR20L160	1,45	2,0	6	160	GREEN	60
CSSR30L130	2,54	3,0	9	130	BLUE	60
CSSR30L160	2,54	3,0	9	160	BLUE	60
CSSR40L130	3,0	4,0	12	130	WHITE	60
CSSR40L160	3,0	4,0	12	160	WHITE	60

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:



SIM MICRO CANNULA “Soft Touch”

Steel suction cannula with suction control handle and atraumatic silicone tip

Revision date: 28.08.2026

Atraumatic silicone tip
made of “Medical Grade” silicone
Shore 60

Why?

- Because reusable suction cannulae with fine lumens are very difficult to clean
- Because disposable suction cannulae are recommended for very high risk procedures, especially neuro and ophthalmic

Characteristics

- Complete line and specials can be easily produced following your indications
- Atraumatic silicone tip, made of “Medical Grade” silicone Shore 60
- Non blocking design
- Cannulae are made of Steel AISI 304
- Single angle
- Internal diameter is uniform over the entire length of the cannula
- Suction control handle made of SBC plastic material
- Matches all commonly used connection tubes
- Colour Code for prompt recognition
- Latex free
- Sterilised by gamma radiation. Single sterile blister packaging



CLASS IIa

0425

RULE 6

STERILE R



Code	Internal Diameter mm.	External Diameter mm.	CH	Cannula length mm.	Handle colour	Pcs. per BOX	Applications
CSSR20X80S	1,45	2,0	6	80	GREEN	60	ENT, Plastics, Otological Surgery & outpatients
CSSR20X130S	1,45	2,0	6	130	GREEN	60	General, Plastics, Obs & Gynae, Paediatrics, Neuro, Oral & Max Fax
CSSR20X160S	1,45	2,0	6	160	GREEN	60	General, Plastics, Obs & Gynae, Paediatrics, Neuro, Oral & Max Fax
CSSR30X80S	2,54	3,0	9	80	BLUE	60	ENT, Plastics
CSSR30X130S	2,54	3,0	9	130	BLUE	60	General, Plastics, Obs & Gynae, Paediatrics, Neuro, Oral & Max Fax
CSSR30X160S	2,54	3,0	9	160	BLUE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics
CSSR40X80S	3,0	4,0	12	80	WHITE	60	ENT, Plastics
CSSR40X130S	3,0	4,0	12	130	WHITE	60	General, Plastics, Obs & Gynae, Paediatrics, Neuro, Oral & Max Fax
CSSR40X160S	3,0	4,0	12	160	WHITE	60	General, Plastics, Paediatrics, Neuro, Oral & Max Fax, fine Orthopaedics

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu


E-Mail: info@medeurope.eu

Distributed by:



SIM MICRO HANDLE

Aspiration control handle for needles / cannulae



- Better hold and manoeuvrability
- Complete with aspiration control
- Made of plastic material
- Latex free
- Single sterile packaging
- Sterilised by gamma radiation
- Non invasive

REF. **SMH**

Box 40 Pcs.



CLASS Is

RULE 1

0425

STERILE R



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simitalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:



SIM MICRO CANNULA for microlaryngoscopy


Steel suction cannula with suction control handle

- Aspiration control hole has a diameter of 1,5 mm., allowing the surgeon a slower aspiration response.
- Whenever the hole is opened the aspiration on the tip will slowly decrease. Even with a very high aspiration pressure the cannula makes only a low amount of noise.
- Particularly recommended for DIRECT SUSPENSION MICROLARYNGOSCOPY
- Complete line and specials can be easily produced following your indications
- Rounded edges, atraumatic distal end of lumen
- Non blocking design
- Cannulae are made of Steel AISI 304
- Single angle
- Internal diameter is uniform over the entire length of the cannula
- Suction control handle made of SBC plastic material
- Matches all commonly used connection tubes
- Colour Code for prompt recognition
- Latex free
- Sterilised by gamma radiation
- Single sterile packaging

CE
0425

CLASS IIa
RULE 6

STERILE R



Code	Internal Diameter mm.	External Diameter mm.	CH	Cannula length mm.	Handle colour	Pcs. per BOX
CSSR30X330	2,54	3,0	9	330	BLUE	30

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu



Distributed by:



SIM MICRO CANNULA “Slow Reaction”

Steel suction cannula with drop-shaped suction control handle


Revision date: 28.08.2026



- By its drop shaped suction interrupter, FUKUSHIMA-Design suction cannula ensures excellent fine control of the suction power
- Rounded edges, atraumatic distal end of lumen
- Cannulae are made of Steel AISI 304
- Internal diameter is uniform over the entire length of the cannula
- Suction control handle made of SBC plastic material
- Matches all commonly used connection tubes
- Colour code for prompt recognition
- Latex free
- Sterilised by gamma radiation
- Packaging: outer sterile blister and inner protective bag

Code	Internal Diameter mm.	External Diameter mm.	CH	Cannula length mm.	Handle colour	Pcs. per BOX
CSSR14G090*	0,96	1,4	4	90	WHITE	60
CSSR14G130*	0,96	1,4	4	130	WHITE	60
CSSR20G80*	1,45	2,0	6	80	GREEN	60
CSSR20G130*	1,45	2,0	6	130	GREEN	60
CSSR20G160*	1,45	2,0	6	160	GREEN	60
CSSR30G80*	2,54	3,0	9	80	BLUE	60
CSSR30G130*	2,54	3,0	9	130	BLUE	60
CSSR30G160*	2,54	3,0	9	160	BLUE	60
CSSR40G80*	3,0	4,0	12	80	WHITE	60
CSSR40G130*	3,0	4,0	12	130	WHITE	60
CSSR40G160*	3,0	4,0	12	160	WHITE	60

*Product not in stock, available by specific request



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:

