


Retractor for NEUROSURGERY

Hook for Galea



Revision date: 26.01.2017

Description

Hooked retractor for galea for use in neurosurgery.

Characteristics

- Made of AISI 316L steel
- Diameter 1,6mm.
- “Medical grade” silicone elastic attachment
- Latex free device
- Packaged in double sterile envelope
- Gamma ray sterilization
- Box of 50 pcs.



CLASS IIa

RULE 7

0476

STERILE R



Codes

REF MEDHOOK001

Large size, 110mms.

REF MEDHOOK002


Large size, 110mms., coated with “medical grade” silicone

REF MEDHOOK003

Small size, 100mms.

REF MEDHOOK004

Small size, 100mms., coated with “medical grade” silicone



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:



Retractor for NEUROSURGERY

Double Hook for Galea

Revision date: 22.10.2020



Description



Double hook retractor for galea for use in neurosurgery.

Characteristics

- Made of AISI 301 steel
- Thickness 0,7 mm
- Length 33 mm
- "Medical grade" silicone elastic attachment
- Latex free device
- Packaged in double sterile envelope
- Sterilised by gamma rays
- Box of 50pz.

Code

REF MEDHOOK-005
(double hook)



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy
Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845
www.simitalia.it www.medeurope.eu
E-Mail: info@medeurope.eu

Distributed by:

