



# Med Stomach Tube


Gastric calibration tube



The device is indicated for use in surgical procedures unique to the stomach. Introduced by mouth of the patient it allows, through its asymmetrical distal balloon, the movement of tissues, facilitating the surgeon involved in the intervention in laparoscopic or open surgery. The holes on the distal portion of atraumatic tip favor the complete suction.

## Characteristics

- Material: medical grade silicone
- With radiopaque line
- With or without asymmetrical distal balloon
- Length 800 mm. or 1000mm
- Sterilised by gamma rays
- Latex and Phthalates free
- Double packed



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy  
 Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845  
 www.simitalia.it www.medeurope.eu  
 E-Mail: info@medeurope.eu

Distributed by:

**arfi**

<b>CODICE / REF</b>	<b>CH</b>	<b>Shore</b>	<b>Length</b>	<b>Modello</b>
MST-34-BL-01-080*	34	80	80 cm	with asymmetrical distal balloon
MST-36-BL-01-080	36	80	80 cm	with asymmetrical distal balloon
MST-38-BL-01-080	38	80	80 cm	with asymmetrical distal balloon
MST-40-BL-01-080*	40	80	80 cm	with asymmetrical distal balloon
MST-42-BL-01-080*	42	80	80 cm	with asymmetrical distal balloon
MST-44-BL-01-080*	44	80	80 cm	with asymmetrical distal balloon
MST-34-AF-01-080*	34	80	80 cm	-
MST-36-AF-01-080*	36	80	80 cm	-
MST-38-AF-01-080	38	80	80 cm	-
MST-40-AF-01-080*	40	80	80 cm	-
MST-42-AF-01-080*	42	80	80 cm	-
MST-44-AF-01-080*	44	80	80 cm	-
MST-34-LL-01-080*	34	57	80 cm	-
MST-36-LL-01-080*	36	57	80 cm	-
MST-38-LL-01-080*	38	57	80 cm	-
MST-40-LL-01-080	40	57	80 cm	-
MST-42-LL-01-080*	42	57	80 cm	-
MST-44-LL-01-080*	44	57	80 cm	-
MST-34-LL-02-100*	34	57	100 cm	with black tip
MST-36-LL-02-100*	36	57	100 cm	with black tip
MST-38-LL-02-100*	38	57	100 cm	with black tip
MST-40-LL-02-100*	40	57	100 cm	with black tip
MST-42-LL-02-100*	42	57	100 cm	with black tip
MST-44-LL-02-100*	44	57	100 cm	with black tip

\* Product not in stock; available by specific request

\* Product not in stock, available by specific request

Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy  
 Tel. +39 051.817902 +39 051.817850 Fax +39 051.817845  
 www.simitalia.it www.medeurope.eu  
 E-Mail: info@medeurope.eu

Distributed by:

