




ELIZABETH

Polyethylene expander for tissue traction or dilation during surgical procedures

Revision date: 31.08.2026



Description

The main component of the medical device, necessary to maintain tissue expansion. It features a studded surface to improve anchoring and stability (Stability System).

Characteristics

- Made of HDPE resin.
- Latex free.
- Packaged in internal blister and external pouch.
- Sterilized by gamma rays.
- Box of 5 pcs.



Codes

REF SMA-SYST-S

Small studded expander

REF SMA-SYST-M

Medium studded expander

REF SMA-SYST-L

Large studded expander



0425

CLASS IIa

RULE 7

STERILE R



Manufactured by:



Via Gramsci, 9 - 40018 San Pietro in Casale (BO) - Italy

Tel. +39 051.817902 +39 051.817850

www.simalia.it www.medeurope.eu

E-Mail: info@medeurope.eu

Distributed by:



DISPOSITIVI MEDICALI DAL 1972